

ENGINEERED SHARPS INJURY
PREVENTION

ERGONOMIC WORKPLACE
SAFETY

MEDICATION AND SPECIMEN
LABELING

WRONG-SITE SURGERY
PREVENTION

ANSELL — THE COMPANY YOU TRUST FOR
HAND PROTECTION NOW OFFERS YOU THE
FULL LINE OF SANDEL SAFETY PRODUCTS



SANDEL

**ANTICIPATE THE
UNEXPECTED IN
PERIOPERATIVE SAFETY**

Ansell



HEALTHCARE SAFETY DEVICES THAT ENHANCE PROTECTION FOR HEALTHCARE WORKERS AND PATIENTS

TOGETHER, WE'RE RAISING THE STANDARDS OF PERIOPERATIVE SAFETY

Ansell, a global leader in protection solutions, is proud to extend its reach into the perioperative environment with the complete SANDEL line of healthcare safety devices.

Ansell is dedicated to providing safety solutions that eliminate, prevent, or significantly reduce workplace accidents and resulting injuries before, during, and after surgery.

By collaborating with nurses, surgical teams, and healthcare experts, we have been able to identify important safety issues and determine how the standard of care can be elevated to reduce safety hazards for both patients and healthcare staff. The resulting SANDEL products have been developed to resolve the specific safety and ergonomic needs identified in a cost-effective manner. It is our hope that these innovative healthcare safety devices will deliver the added confidence and peace of mind all healthcare professionals need to perform at their best.

TABLE OF CONTENTS

1. ENGINEERED SHARPS INJURY PREVENTION	3-6
Change-A-Blade®	3
Weighted and Non-Weighted Safety Scalpels	4
Scalpel Disarmer™	4
Neutral Zone	5
Hands-Free Transfer Trays	5
Ordering Information	6
2. ERGONOMIC WORKPLACE SAFETY	7-10
Patient Transfer System	7
Cord Cover	7
Absorbent Floor Pad	8
Step Stools	8
Anti-Fatigue Mats	8
Sit-Stand Stool	8
Ordering Information	9-10
3. MEDICATION AND SPECIMEN LABELING	11-14
Correct Medication Labeling System™	11
Correct Specimen Handling System™	12
Ordering Information	13-14
4. WRONG-SITE SURGERY PREVENTION	15-18
TIME OUT® Hood	15
TIME OUT® Beacon	15
TIME OUT® Sleeve	15
TIME OUT® Markers	16
Correct-Site Surgery Stickers & Tattoos	16
Ordering Information	17-18
5. UTILITY PRODUCTS	19-20
Permanent Markers	19
Cord & Tube Holder	19
Suture Bag	19
Writing Pen with Light	20
Ordering Information	20
6. SAFETY KITS	21-22
TIME OUT® Safety Kits	21
Change-A-Blade® Safety Kits	21
Ordering Information	22



HELPING ELIMINATE THE THREAT OF SHARPS INJURIES

In the United States, 600,000 to 800,000 sharps injuries occur annually.¹ Ansell has anticipated the need for safer ways to handle sharps and responded with the full line of SANDEL Engineered Sharps Injury Prevention (ESIP) products. These devices are designed to eliminate the threat of sharps injuries, potentially reduce the transmission of blood-borne pathogens, and help healthcare workers comply with OSHA, AORN, and ACS.



CHANGE-A-BLADE®

The SANDEL Change-A-Blade® is a disposable scalpel handle with a built-in protective safety shield.

- Compatible with any manufacturer's stainless steel or carbon steel blades
- Weighted to the feel of a conventional handle
- Does not require any change in the surgeon's technique
- Allows for multiple blade changes in the same procedure
- Available in a range of handle sizes compatible with multiple standard blade sizes
- Comes with a removable TIME OUT® Sleeve*



WEIGHTED AND NON-WEIGHTED SAFETY SCALPELS

The SANDEL disposable Weighted and Non-Weighted Safety Scalpels provide a sliding shield to cover the surgical blade and help prevent injuries.

- Available in a wide variety of surgical blade and handle sizes
- Mimic the feel of a standard scalpel
- Color-coded protective safety shields for proper blade identification
- Safety-shield button easily locks and unlocks with either left- or right-handed use
- Come with a removable TIME OUT® Sleeve*



SCALPEL DISARMER™

The SANDEL Scalpel Disarmer™ is a blade and needle counting and disposal box designed to help safely remove and dispose of scalpel blades and needles.

- Two separate, removable halves with adhesive tape to provide safety and stability on any surface
- Easily locks for safe disposal
- Specially designed mechanism for safely removing blades from the Change-A-Blade® or any standard scalpel handle
- Available with a magnet or foam surface for counting needles

*Please see pages 15-16 for more information about TIME OUT®.



ORDERING INFORMATION



NEUTRAL ZONE

The SANDEL Z-Friction Drape™ provides a flexible, non-skid neutral zone on the sterile field for the placement of sharps and minimization of hand-to-hand transfer.

- Distinctive bright orange color
- Adaptable design, allowing it to be cut and placed in two locations
- Non-magnetic, non-slip design holds instruments in place and keeps them from sliding off the sterile field



HANDS-FREE TRANSFER TRAYS

The SANDEL Z-Tray™ and Stretch-A-Tray™ protect the surgical team by providing a designated area for placement and retrieval of sharps, minimizing the risk of sharps injuries.

- Distinctive bright orange color
- Two different designs to accommodate a wide variety of sharps, handles, and instruments
- Includes adhesive strip for use as a non-movable neutral zone where sharps can be safely placed


Product Type	Code Number	Product Name	Description	Packaging	Case Qty.
Change-A-Blade®	2200	Change-A-Blade® #3 Handle	Weighted Disposable Safety Handle (#3), No Blade, (fits any #10, #11, or #15 blades) with TIME OUT® Sleeve (S)	12/Bx 4 Bx/Cs	48
	2200-L	Change-A-Blade® #3L Handle	#3L Weighted Disposable Safety Handle (8.75"), No Blade, (fits any #10, #11, or #15 blades) with TIME OUT® Sleeve (S)	12/Bx 2 Bx/Cs	24
	2204	Change-A-Blade® #4 Handle	Weighted Disposable Safety Handle (#4), No Blade, (fits any #20, #22, #23, or #24) with TIME OUT® Sleeve (S)	12/Bx 4 Bx/Cs	48
Weighted Safety Scalpel	2210-N	Weighted Safety Scalpel	#3 Handle with #10 Blade, Orange Shield with TIME OUT® Sleeve (S)	12/Bx 8 Bx/Cs	96
	2210-L	Weighted Safety Scalpel	#3L Handle (8.75") with #10 Blade, Orange Shield with TIME OUT® Sleeve (S)	12/Bx 2 Bx/Cs	24
	2211-N	Weighted Safety Scalpel	#3 Handle with #11 Blade, Blue Shield with TIME OUT® Sleeve (S)	12/Bx 8 Bx/Cs	96
	2211-L	Weighted Safety Scalpel	#3L Handle (8.75") with #11 Blade, Blue Shield with TIME OUT® Sleeve (S)	12/Bx 2 Bx/Cs	24
	2215-N	Weighted Safety Scalpel	#3 Handle with #15 Blade, Yellow Shield with TIME OUT® Sleeve (S)	12/Bx 8 Bx/Cs	96
	2215-L	Weighted Safety Scalpel	#3L Handle (8.75") with #15 Blade, Yellow Shield with TIME OUT® Sleeve (S)	12/Bx 2 Bx/Cs	24
Non-Weighted Safety Scalpel	2220-N	Weighted Safety Scalpel	#4 Handle with #20 Blade, Orange Shield with TIME OUT® Sleeve (S)	12/Bx 4 Bx/Cs	48
	2230-N	Non-Weighted Safety Scalpel	#3 Handle with #10 Blade, Orange Shield with TIME OUT® Sleeve (S)	12/Bx 8 Bx/Cs	96
	2231-N	Non-Weighted Safety Scalpel	#3 Handle with #11 Blade, Blue Shield with TIME OUT® Sleeve (S)	12/Bx 8 Bx/Cs	96
Scalpel Disarmer™	2241	Scalpel Disarmer™	Scalpel Disarmer™ Magnet & Foam 40 Count (S)	12/Bx 4 Bx/Cs	48
	2242	Scalpel Disarmer™	Scalpel Disarmer™ Magnet & Magnet 30 Count (S)	12/Bx 4 Bx/Cs	48
Neutral Zone	1151	Z-Friction Drape™	Instrument Drape (10"x 16") (S)	40/Bx 2 Bx/Cs	80
Hands-Free Transfer Trays	1145	Z-Tray™	Hands-Free Transfer Tray with Adhesive Strip (S)	12/Bx 2 Bx/Cs	24
	1147	Stretch-A-Tray™	Extendable Hands-Free Transfer Tray with Adhesive Strip (S)	10/Bx 2 Bx/Cs	20

S=sterile
NS=non-sterile



PROTECTING HEALTHCARE STAFF FROM COSTLY WORKPLACE INJURIES

The incidence rate of lost-workday injuries from slips, trips, and falls for US hospital employees is 90% greater than the average rate for all other private industries *combined*. These types of injuries are the second most common cause of lost-workday injuries in hospitals.² Additionally, more than \$15 to \$20 billion is spent in workers' compensation costs for musculoskeletal disorders each year.³ SANDEL safety solutions offer cost-effective ways to help prevent these common injuries and comply with OSHA regulations for workplace safety.



PATIENT TRANSFER SYSTEM

The SANDEL Z-Slider® Patient Transfer & Repositioning Sheet helps prevent disabling back injuries to staff by decreasing physical strain to the back, shoulders, neck, and arms.

- Supports proper patient handling and ergonomic techniques
- Minimizes the potential for patient skin tears, abrasions, and bruising
- Disposable, reducing the potential for cross-contamination
- Single-patient use, designed for multiple uses with the same patient
- No weight limit



CORD COVER

The SANDEL Trip-No-More™ is a bright orange disposable cord cover designed to keep cords and tubes in place on the floor and reduce trips and falls.

- Adhesive strips keep it firmly in place
- Easy to apply and remove—leaves no residue on the floor
- Distinctive orange color
- Available in a variety of lengths



ABSORBENT FLOOR PAD

The SANDEL DriFloor™ Absorbent Floor Pad collects fluids to prevent slips and falls and improve the efficiency of cleanup, reducing costly turnover time.

- Eliminates the unsanitary and inappropriate use of blankets or towels on the floor
- Fluid-proof backing keeps the floor dry underneath the pad
- Available in a variety of sizes



STEP STOOLS

The SANDEL Ergo-Step™ Stools are lightweight, stackable, "next generation" replacements for traditional, heavy, metal step stools.

- Lightweight (~8lbs) and easy to handle
- Weight capacity of approximately 500 lbs
- Anti-skid feet keep it firmly in place
- Distinctive bright orange color
- Can be interlocked to create a safe work platform



ANTI-FATIGUE MATS

SANDEL Ergo-Stand EZ™ Anti-Fatigue Mats reduce stress and strain on joints caused by standing in one place for an extended period of time.

- Lightweight and reusable
- Available in a variety of sizes
- Beveled on all sides to prevent tripping
- Can be placed on step stools during long procedures
- Disposable cover available for solo mat to keep it clean and reduce turnover time



SIT-STAND STOOL

The Standel™ Sit-Stand Stool helps reduce fatigue and relieve pressure on the feet, legs, and back due to standing in the same position during long procedures.

- Designed for use in multiple locations such as the OR or nurses' stations
- Adjustable height from 20" to 30"
- Weight capacity of 330 lbs



ERGONOMIC WORKPLACE SAFETY



ORDERING INFORMATION

Product Type	Code Number	Product Name	Description	Packaging	Case Qty.
Patient Transfer System	2102	Z-Slider®	Patient Transfer & Repositioning Sheet (on a roll) (17½" x 5" x 5")	30/Bx 2 Bx/Cs	60
	2103	Z-Slider®	Patient Transfer & Repositioning Sheet Large Dispenser Box (11⅞" x 20⅞" x 1⅞")	15/Bx 4 Bx/Cs	60
	2104	Z-Slider®	Patient Transfer & Repositioning Sheet Small Dispenser Box (11" x 9¾" x 2")	10/Bx 6 Bx/Cs	60
	2103-H	Z-Slider® Clear Wall Holder	Point-of-Use Clear Wall Mounting Holder for the 2103 Dispenser Box (12¼" x 20⅞" x 1⅞")	5/Cs	5
	2104-H	Z-Slider® Clear Wall Holder	Point-of-Use Clear Wall Mounting Holder for the 2104 Dispenser Box (11⅞" x 10" x 2⅔")	5/Cs	5
	2103-WB	Z-Slider® Metal Mounting Bracket	Mounting Bracket for use with 2103 (11" x 2")	10/Bx or sold individually	10
Cord Cover	2302	Trip-No-More™ Cord Alert	Cord Alert (8" x 24")	25/Bag 3 Bags/Cs	75
	2303	Trip-No-More™ Cord Alert	Cord Alert (8" x 36")	25/Bag 3 Bags/Cs	75
	2309	Trip-No-More™ Cord Alert	Cord Alert (8" x 125' roll)	1 Bag of 1 Roll/Cs	1
	2312	Trip-No-More™ Cord Alert	Cord Alert (8" x 24") (Glow-in-the-dark)	75/Bag 1 Bag/Cs	75
	2313	Trip-No-More™ Cord Alert	Cord Alert (8" x 36") (Glow-in-the-dark)	75/Bag 1 Bag/Cs	75
Absorbent Floor Pad	2419	DriFloor™ Pad	Absorbent Floor Pad (16½" x 20") with Backing	25/Bag 3 Bags/Cs	75
	2440	DriFloor™ Pad	Absorbent Floor Pad (33" x 40") with Backing	20/Bag 1 Bag/Cs	20
	2449	DriFloor™ Roll	Absorbent Floor Pad on a Roll (33" x 100') with Backing	1 Bag of 1 Roll/Cs	1

Product Type	Code Number	Product Name	Description	Packaging	Case Qty.
Step Stools	1170	Ergo-Step™ Stool	2 Tops and 2 Bottoms	4/Cs	4
	1170-T	Ergo-Step™ Stool	4 Tops	4/Cs	4
	1170-B	Ergo-Step™ Stool	4 Bottoms	4/Cs	4
	1170-C	Ergo-Step™ Stool Cart	Holds a maximum of 6 Step Stools	1/Cs	1
Anti-Fatigue Mats	2317	Solo Ergo-Stand EZ™	Anti-Fatigue Mat—Reusable (13" x 17") (NS)	10/Bx 2 Bx/Cs	20
	2317-C	Solo Ergo-Stand EZ™ Cover	Disposable Cover for 2317 (14" x 18") (NS)	100/Bag 2 Bags/Cs	200
	2317-BC	Solo Ergo-Stand EZ™	Solo Anti-Fatigue Mat—Reusable (13" x 17") (NS)	10/Bag 1 Bag/Cs	10
	2317-S	Solo Ergo-Stand EZ™	Solo Anti-Fatigue Mat—Reusable (13" x 17") (NS)	1/Cs	1
	2317-S3	Solo Ergo-Stand EZ™	Solo Anti-Fatigue Mat—Reusable (13" x 17") (NS)	3/Cs	3
	2324	Ergo-Stand EZ™	Anti-Fatigue Mat—Reusable (24" x 24") (NS)	8/Cs	8
	2348	Ergo-Stand EZ™	Anti-Fatigue Mat—Reusable (24" x 48") (NS)	4/Cs	4
Sit-Stand Stool	2300	Standel™ Sit-Stand	Adjustable Sit-Stand Stool	1/Cs	1

S=sterile
NS=non-sterile



ENHANCING COMMUNICATION AND REDUCING ERRORS

39% of medication errors are due to naming, labeling, or packaging.⁴ Medication errors cause at least one death every day in the United States and injure approximately 1.3 million people annually.⁵ The SANDEL Correct Medication Labeling System™ provides the tools you need to accurately label medications to help prevent potentially costly errors and keep your facility in compliance with The Joint Commission's requirements.



CORRECT MEDICATION LABELING SYSTEM™

Ansell provides a wide range of medication labeling kits customized specifically for each environment—from ophthalmic to labor & delivery—designed for clarity and functionality.

- Kits can include a variety of medication labels, markers, flags, medication cups, and more to suit your needs
- Proprietary flag system increases the visibility of labels and provides a color-coded strip for the syringe
- Labels are moisture- and smear-resistant and will not tear or fall off when submerged in fluid
- Labels are tailored for individual specialties and preprinted with department-specific medication names
- Kits can be customized based on your department's specific needs

IMPROVING ACCURACY AND PATIENT CARE

Hospitals experience nearly one error per 1,000 specimen labels.⁶ Nearly 30% of these errors are mislabeled specimens, the most common error type. This is especially problematic since these types of mistakes are the ones that often have the greatest impact on patient safety and healthcare costs.⁶ The SANDEL Correct Specimen Handling System™ is designed to improve clarity and help eliminate these errors to improve overall accuracy and safety.



CORRECT SPECIMEN HANDLING SYSTEM™

The SANDEL Correct Specimen Handling System™ provides a multitude of options to label and organize specimens.

- Labels allow for clear, lasting identification of specimens to help reduce errors
- Distinctive bright orange color
- Containers available in both sterile and non-sterile options



MEDICATION AND SPECIMEN LABELING



ORDERING INFORMATION

Product Type	Code Number	Product Name	Description	Packaging	Case Qty.
Correct Medication Labeling System™—OR	3302	Correct Medication Labeling System™ Designed for the OR	(4) Flags, (1) 4-in-1 Marker™, (2) Sheets of 24 Preprinted Labels, (2) Med Cups, (1) Specimen Zone™ Strip (4 count) (S)	25/Bx 2 Bx/Cs	50
	3303-P	Correct Medication Labeling System™ Designed for the OR	(4) Flags, (1) 4-in-1 Marker™, (2) Sheets of 24 Preprinted Labels, (1) Specimen Zone™ Strip (4 count) (S)	50/Bx 2 Bx/Cs	100
	3304	Correct Medication Labeling System™ Designed for the OR	(4) Flags, (2) Sheets of 24 Preprinted Labels, (1) Specimen Zone™ Strip (4 count) (S)	50/Bx 1 Bx/Cs	50
	3310-P	Correct Medication Labeling System™ Designed for the OR	(2) Flags, (1) 4-in-1 Marker™, (2) Sheets of 24 Preprinted Labels, (1) Specimen Zone™ Strip (4 count) (S)	50/Bx 2 Bx/Cs	100
	3311	Correct Medication Labeling System™ Designed for the OR	(2) Flags, (1) Petite Permanent Ink Marker, (2) Sheets of 16 Preprinted Labels (S)	100/Bx 1 Bx/Cs	100
	3312	Correct Medication Labeling System™ Designed for the OR	(2) Flags, (2) Sheets of 24 Preprinted Labels, (1) Specimen Zone™ Strip (4 count) (S)	50/Bx 1 Bx/Cs	50
	3330-P	Correct Medication Labeling System™ Designed for the OR	(2) Sheets of 24 Preprinted Labels (S)	100/Bx 2 Bx/Cs	200
	3331-P	Correct Medication Labeling System™ Designed for the OR	(2) Sheets of 16 Preprinted Labels (S)	100/Bx 2 Bx/Cs	200
	3333	Correct Medication Labeling System™ Designed for the OR	(1) Petite Permanent Ink Marker, (2) Sheets of 16 Preprinted Labels (S)	100/Bx 1 Bx/Cs	100
	3334-P	Correct Medication Labeling System™ Designed for the OR	(1) 4-in-1 Marker™, (2) Sheets of 24 Preprinted Labels, (1) Specimen Zone™ Strip (4 count) (S)	50/Bx 2 Bx/Cs	100
Correct Medication Labeling System™—General Use	3300-1	Correct Medication Labeling System™ Designed for General Use	(2) Flags, (1) 4-in-1 Marker™, (1) Sheet of 16 Preprinted Labels with Prompts for Medication Name and Strength (S)	50/Bx 1 Bx/Cs	50
	3338-T	Correct Medication Labeling System™ Designed for General Use	(1) Sheet of 8 Preprinted Labels with Prompts. Compliant with New 2010 Joint Commission Requirements (S)	100/Bx 2 Bx/Cs	200
Correct Medication Labeling System™—Surgery Centers	3300	Correct Medication Labeling System™ Designed for Surgery Centers	(2) Flags, (1) 4-in-1 Marker™, (1) Sheet of 16 Preprinted Labels, (1) Specimen Zone™ Strip (4 count) (S)	50/Bx 1 Bx/Cs	50
	3300-L	Correct Medication Labeling System™ Designed for Surgery Centers	(1) Sheet of 16 Preprinted Labels (S)	100/Bx 1 Bx/Cs	100
Correct Medication Labeling System™—Cath Labs	3301	Correct Medication Labeling System™ Designed for Cath Labs	(4) Flags, (1) 4-in-1 Marker™, (2) Sheets of 16 Preprinted Labels (S)	50/Bx 1 Bx/Cs	50
	3301-L	Correct Medication Labeling System™ Designed for Cath Labs	(2) Sheets of 16 Preprinted Labels (S)	100/Bx 1 Bx/Cs	100

Product Type	Code Number	Product Name	Description	Packaging	Case Qty.
Correct Medication Labeling System™—Interventional Radiology	3350	Correct Medication Labeling System™ Designed for Interventional Radiology	(2) Flags, (1) 4-in-1 Marker™, (2) Sheets of 16 Preprinted Labels (S)	50/Bx 1 Bx/Cs	50
	3350-L	Correct Medication Labeling System™ Designed for Interventional Radiology	(2) Sheets of 16 Preprinted Labels (S)	100/Bx 1 Bx/Cs	100
Correct Medication Labeling System™—Ophthalmic Surgery	3360	Correct Medication Labeling System™ Designed for Ophthalmic Surgery	(1) Flag, (1) Single-Tip "No Roll" Permanent Ink Marker, (1) Sheet of 16 Preprinted Labels (S)	50/Bx 1 Bx/Cs	50
	3360-L	Correct Medication Labeling System™ Designed for Ophthalmic Surgery	(1) Sheet of 16 Preprinted Labels (S)	100/Bx 1 Bx/Cs	100
Correct Medication Labeling System™—Pain Management	3370	Correct Medication Labeling System™ Designed for Pain Management	(2) Flags, (1) 4-in-1 Marker™, (2) Sheets of 16 Preprinted Labels (S)	50/Bx 1 Bx/Cs	50
	3370-L	Correct Medication Labeling System™ Designed for Pain Management	(2) Sheets of 16 Preprinted Labels (S)	100/Bx 1 Bx/Cs	100
Correct Medication Labeling System™—Labor & Delivery	3390	Correct Medication Labeling System™ Designed for Labor & Delivery	(2) Flags, (1) 4-in-1 Marker™, (1) Sheet of 16 Preprinted Labels (S)	50/Bx 1 Bx/Cs	50
	3390-L	Correct Medication Labeling System™ Designed for Labor & Delivery	(2) Sheets of 16 Preprinted Labels (S)	100/Bx 1 Bx/Cs	100
ID Flags	3305	Flags with Matching Syringe ID Strip	4-Flag Pack (1 ea. Orange, Red, Yellow, Blue) (S)	50/Bx 1 Bx/Cs	50
Correct Specimen Handling System™	3450	Correct Specimen Zone™	Orange Specimen Zone (4 Count), (1) Petite Permanent Marker (S)	50/Bx 2 Bx/Cs	100
	3455	Orange Specimen Labels	Orange Specimen Labels (2" x 4") on Roll (NS)	500/Roll 1 Roll/Cs	500
	3460	Orange Specimen Wrap	Orange Specimen Wrap—Cotton (23½" x 16½") (S)	24/Bx 4 Bx/Cs	96
	3460-NS	Orange Specimen Wrap	Orange Specimen Wrap—Cotton (23½" x 16½") (NS)	100/Bag 1 Bag/Cs	100
	3470	Orange Cup with Lid	120cc Cup with Orange Lid (S)	25/Bx 2 Bx/Cs	50
	3475	Shaped Orange Specimen Bowl	Shaped Orange Specimen Bowl (S)	20/Bag 1 Bag/Cs	20
	3475-NS	Shaped Orange Specimen Bowl	Shaped Orange Specimen Bowl (NS)	100/Bag 1 Bag/Cs	100
	3480	Orange Specimen Bag	Large Orange Specimen Bag (33" x 39") (S)	25/Bag 2 Bags/Cs	50

S=sterile
NS=non-sterile



HELPING TO REDUCE SURGICAL ERRORS

The Joint Commission's Universal Protocol, as part of their National Patient Safety Goals, requires taking a "time out" to verify the site, procedure, and patient prior to the start of surgery—helping to avoid potential mistakes.⁷ SANDEL TIME OUT® products serve as prominent reminders to the surgical team to comply with these requirements and help reduce the incidence of errors.



TIME OUT® HOOD

The TIME OUT® Hood can be placed over instruments, to be removed before use, as a reminder to take the "time out" before beginning a procedure.

- Distinctive bright orange color
- Reusable hood can be placed over any scalpel or other instrument
- Can be sterilized in the instrument tray



TIME OUT® BEACON

The TIME OUT® Beacon can be placed over the surgical site, Mayo stand, or back table as a barrier that must be removed before the incision is made to remind the surgical team to take a "time out."

- Distinctive bright orange color
- Available in cotton or a disposable, non-woven biodegradable material
- Can be sterilized with steam or ETO



TIME OUT® SLEEVE

Slip the TIME OUT® Sleeve over a standard scalpel handle before setup or place it on the Mayo stand or back table as a prompt for the surgical team to take the "time out."

- Distinctive bright orange color
- Included with SANDEL Safety Scalpels, Change-A-Blade®, and 4-in-1 Marker™
- Comes with a ruler on the side panel

HELPING TO ENSURE ACCURACY

It is estimated that 1,500 to 2,500 wrong-site incidents occur each year in the United States alone.⁷ SANDEL Wrong-Site Surgery Prevention products provide a variety of ways to mark the correct surgical site and help eliminate mistakes.



TIME OUT® MARKERS

SANDEL TIME OUT® Skin Markers are designed to mark the correct surgical site on a patient's skin.

- Distinctive bright orange color
- Triangular barrel prevents rolling
- Available in both single- and dual-tip with multiple thickness options



CORRECT-SITE SURGERY STICKERS & TATTOOS

Correct-Site Surgery Stickers & Tattoos can be placed on the patient to mark the correct surgical location.

- Green and white coloring stands out against all skin colors
- Available in multiple sizes



WRONG-SITE SURGERY PREVENTION



ORDERING INFORMATION

Product Type	Code Number	Product Name	Description	Packaging	Case Qty.
TIME OUT® Hood	711	TIME OUT® Hood	TIME OUT® Hood—Reusable (NS)	10/Bag 1 Bag/Cs	10
TIME OUT® Beacon	811	TIME OUT® Beacon	TIME OUT® Reminder—Disposable (Cotton) (S)	24/Bx 4 Bx/Cs	96
	811-E	TIME OUT® Beacon	TIME OUT® Reminder—Disposable (Non-Woven) Biodegradable (S)	100/Bx 2 Bx/Cs	200
	811-NS	TIME OUT® Beacon	TIME OUT® Reminder—Disposable (Cotton) (NS)	100/Bag 2 Bags/Cs	200
	811-ENS	TIME OUT® Beacon	TIME OUT® Reminder—Disposable (Non-Woven) Biodegradable (NS)	250/Bag 4 Bags/Cs	1000
TIME OUT® Sleeve	911	TIME OUT® Sleeve	TIME OUT® Sleeve—Disposable (S)	300/Bx 1 Bx/Cs	300
	911-NS	TIME OUT® Sleeve	TIME OUT® Sleeve—Disposable (NS)	500/Bag 1 Bag/Cs	500

Product Type	Code Number	Product Name	Description	Packaging	Case Qty.	
TIME OUT® Markers	1011-NNS	SANDEL Petite Skin Marker™	Petite Single Tip “No Roll” Skin Marker (4” in length) (NS)	50/Bag 10 Bags/Cs	500	
	1019	SANDEL Skin Marker™	Single Tip “No Roll” Skin Marker (S)	25/Bx 4 Bx/Cs	100	
	1019-NNS	SANDEL Skin Marker™	Single Tip “No Roll” Skin Marker (NS)	200/Bx 2 Bx/Cs	400	
	1019-R	SANDEL Skin Marker™ with flexible ruler	Single Tip “No Roll” Skin Marker (S) with 6” flexible ruler	25/Bx 4 Bx/Cs	100	
	1031	SANDEL 2-in-1 Marker™	Dual Ink (Permanent & Skin Marker Ink) Regular Tip Marker with TIME OUT® Sleeve & Ruler (S)	25/Bx 4 Bx/Cs	100	
	1031-L	SANDEL 2-in-1 Marker™	Dual Tip/Dual Ink (Permanent & Skin Marker Ink) Marker with 8 Labels & Ruler (S)	25/Bx 4 Bx/Cs	100	
	1031-NNS	SANDEL 2-in-1 Marker™	Dual Ink (Permanent & Skin Marker Ink) Regular Tip Marker with TIME OUT® Sleeve & Ruler (NS)	200/Bx 2 Bx/Cs	400	
	1041	SANDEL 4-in-1 Marker™	Dual Tip/Dual Ink (Permanent & Skin Marker Ink) Chiseled Tip Marker with TIME OUT® Sleeve & Ruler (S)	25/Bx 4 Bx/Cs	100	
	1041-L	SANDEL 4-in-1 Marker™	Dual Tip/Dual Ink (Permanent & Skin Marker Ink) Chiseled Tip Marker with 8 Labels, TIME OUT® Sleeve & Ruler (S)	25/Bx 4 Bx/Cs	100	
	1041-NNS	SANDEL 4-in-1 Marker™	Dual Tip/Dual Ink (Permanent & Skin Marker Ink) Chiseled Tip Marker with TIME OUT® Sleeve & Ruler (NS)	100/Bag 4 Bags/Cs	400	
	1041-LR	SANDEL 4-in-1 Marker™ with labels & flexible ruler	Dual Tip/Dual Ink (Permanent & Skin Marker Ink) Marker with TIME OUT Sleeve & 6” Flexible Ruler & 8 Labels (S)	25/Bx 4 Bx/Cs	100	
	Correct-Site Surgery Stickers & Tattoos	3175	Correct-Site™ Stickers	Large (3” diameter) Green & White Circle with Signature Line for Patient (NS)	100/Bag 1 Bag/Cs	100
		3176	Correct-Site™ Stickers	Small (1” diameter) Green & White Circle (NS)	100/Bag 1 Bag/Cs	100
3177		Correct-Site™ Wristband	Green Wristband (NS)	100/Bag 1 Bag/Cs	100	
3186		Correct-Site™ Tattoo	Small (1” diameter) Green & White Circle (NS)	50/Bag 5 Bags/Cs	250	

S=sterile
NS=non-sterile



KEEP THE STERILE FIELD ORGANIZED AND EFFICIENT

The SANDEL line offers a cohesive set of products that help maximize organization and safety in the surgical environment.



PERMANENT MARKERS

SANDEL Permanent Markers with a “no-roll” barrel can be used for clear labeling of medications.

- Distinctive bright orange color
- Triangular barrel prevents rolling
- Available in both sterile and non-sterile versions
- Can be purchased individually or as part of a SANDEL Correct Medication Labeling System™



CORD & TUBE HOLDER

The SANDEL Cord & Tube Holder is designed to keep cords and tubes organized on the sterile field.

- Distinctive bright orange color
- Adhesive backing keeps it firmly in place



SUTURE BAG

Keep materials organized with the SANDEL Z-Suture Saddle Bag™, a clear plastic bag that allows for easy identification of its contents.

- Can be draped over the back table or the Mayo stand
- Perforated so the bags can be hung in separate locations



WRITING PEN WITH LIGHT

For any healthcare professional who needs to take notes in a dim environment, the Z-Write Lite™ provides a focused light to increase accuracy of documentation without disturbing the surrounding area.

- Distinctive bright orange color
- Available in sterile and non-sterile versions

ORDERING INFORMATION

Product Type	Code Number	Product Name	Description	Packaging	Case Qty.
Permanent Markers	1000-NNS	SANDEL Permanent Marker	Single Tip “No Roll” Permanent Marker, Black Ink (NS)	100/Bag 4 Bags/Cs	400
	1001-L	SANDEL Permanent Marker	Single Tip “No-Roll” Permanent Marker, Black Ink with 8 Labels Preprinted with Prompts for Medication Name & Strength (S)	25/Bx 4 Bx/Cs	100
	1002-NNS	SANDEL Petite Permanent Marker	Petite Single Tip “No-Roll” Permanent Marker, Black Ink (4” in length) (NS)	50/Bag 10 Bags/Cs	500
Cord & Tube Holder	1117	Z-Cord Holder™	Sterile Holder for Tubes and Cords (S)	50/Bx 2 Bx/Cs	100
Suture Bags	1130	Z-Suture Saddle Bag™	Clear Plastic 20” Saddle Bag for use on Mayo Stand or Back Table (S)	50/Bx 2 Bx/Cs	100
	1134	Z-Suture Saddle Bag™	Clear Plastic 24” Saddle Bag for use on Mayo Stand or Back Table (S)	50/Bx 2 Bx/Cs	100
Writing Pen with Light	1050	Z-Write Lite™	Writing Pen with Light (S)	20/Bx 4 Bx/Cs	80
	1050-NS	Z-Write Lite™	Writing Pen with Light (NS)	50/Bx 2 Bx/Cs	100
	1050-BC	Z-Write Lite™ (Sterile)	Writing Pen with Light (S)	20/Bx 1 Bx/Cs	20
	1050-BCNS	Z-Write Lite™	Writing Pen with Light (NS)	25/Bag 1 Bag/Cs	25

S=sterile
NS=non-sterile



CONSOLIDATE ALL YOUR SAFETY PRODUCTS

SANDEL Safety Kits put everything you need in one place, helping to reduce costs, save time, and reduce waste. A wide variety of options are available to fit your unique needs.



TIME OUT® SAFETY KITS*

TIME OUT® Safety Kits contain everything you need to create a safer surgical environment including helpful tools such as the Weighted Safety Scalpel, Z-Friction Drape™, Stretch-A-Tray™, Colored ID Flags, Medication Labels, 4-in-1 Dual Tip Marker, and more.



CHANGE-A-BLADE® SAFETY KITS

Change-A-Blade® Safety Kits combine the shielded scalpel handle with a number of other products including the Scalpel Disarmer™, Colored ID Flags, Medication Labels, 4-in-1 Dual Tip Marker, and more.

ORDERING INFORMATION

Product Type	Code Number	Product Name	Description	Packaging	Case Qty.
Safety Kits	K-911	TIME OUT® Safety Kit	(1) Z-Friction Drape™ (10" x 16"), (1) Stretch-A-Tray™, (1) Specimen Zone™ Strip (4 count), (4) Color Flags, (2) Sheets of 24 OR Labels, (1) 4-in-1 Dual Tip Marker (S)	12/Bx 2 Bx/Cs	24
	K-911-10	TIME OUT® Safety Kit	(1) #10 Weighted Safety Scalpel, (1) Z-Friction Drape™ (10" x 16"), (1) Z-Tray™, (1) Specimen Zone™ Strip (4 count), (4) Color Flags, (2) Sheets of 24 OR Labels, (1) 4-in-1 Dual Tip Marker (S)	12/Bx 2 Bx/Cs	24
	K-911-11	TIME OUT® Safety Kit	(1) #11 Weighted Safety Scalpel, (1) Z-Friction Drape™ (10" x 16"), (1) Z-Tray™, (1) Specimen Zone™ Strip (4 count), (4) Color Flags, (2) Sheets of 24 OR Labels, (1) 4-in-1 Dual Tip Marker (S)	12/Bx 2 Bx/Cs	24
	K-911-15	TIME OUT® Safety Kit	(1) #15 Weighted Safety Scalpel, (1) Z-Friction Drape™ (10" x 16"), (1) Z-Tray™, (1) Specimen Zone™ Strip (4 count), (4) Color Flags, (2) Sheets of 24 OR Labels, (1) 4-in-1 Dual Tip Marker (S)	12/Bx 2 Bx/Cs	24
	2251	Change-A-Blade® Safety Kit	(1) Change-A-Blade®, (1) Scalpel Disarmer™ with Magnet & Foam (S)	12/Bx 4 Bx/Cs	48
	2252	Change-A-Blade® Safety Kit	(1) Change-A-Blade®, (1) Scalpel Disarmer™ with Magnet & Magnet (S)	12/Bx 4 Bx/Cs	48
	2261	Change-A-Blade® Safety Kit	(1) Change-A-Blade®, (1) Scalpel Disarmer™ with Magnet & Foam, (2) Sheets of 16 OR Labels, (1) 4-in-1 Dual Tip Marker (S)	12/Bx 4 Bx/Cs	48
	2262	Change-A-Blade® Safety Kit	(1) Change-A-Blade®, (1) Scalpel Disarmer™ with Magnet & Foam, (2) Sheets of 16 OR Labels, (1) 4-in-1 Dual Tip Marker (S)	12/Bx 4 Bx/Cs	48
	2271	Change-A-Blade® Safety Kit	(1) Change-A-Blade®, (1) Scalpel Disarmer™ with Magnet & Foam, (2) Sheets of 16 OR Labels, (1) 4-in-1 Dual Tip Marker, (4) Color Flags (S)	12/Bx 4 Bx/Cs	48
	2271-B	Change-A-Blade® Safety Kit	(1) Change-A-Blade®, (1) Scalpel Disarmer™ with Magnet & Foam, (2) Sheets of 16 OR Labels, (1) 4-in-1 Dual Tip Marker, (4) Color Flags, (1) TIME OUT® Beacon (S)	12/Bx 4 Bx/Cs	48
	2272	Change-A-Blade® Safety Kit	(1) Change-A-Blade®, (1) Scalpel Disarmer™ with Magnet & Magnet, (2) Sheets of 16 OR Labels, (1) 4-in-1 Dual Tip Marker, (4) Color Flags (S)	12/Bx 4 Bx/Cs	48

S=sterile
NS=non-sterile

*Please see pages 15-16 for more information about TIME OUT®.



ANTICIPATE THE UNEXPECTED IN PERIOPERATIVE SAFETY

ENGINEERED SHARPS INJURY PREVENTION

ERGONOMIC WORKPLACE SAFETY

MEDICATION AND SPECIMEN LABELING

WRONG-SITE SURGERY PREVENTION

To place your order, visit www.sandelmedical.com or contact your regional office:

Ansell Sandel Medical Solutions LLC
19736 Dearborn Street
Chatsworth, CA 91311
USA
Phone (818) 534-2500
Toll Free (866) 764-3327
Fax (818) 534-2510
E-mail: sandel@ansell.com

To obtain information on the Ansell full line of medical gloves, visit www.ansell.com or contact your regional office:

Ansell Healthcare
200 Schulz Drive
Red Bank, NJ 07701
Phone (732) 345-5400
Toll Free (800) 952-9916
Fax (732) 219-5114
E-mail: info@ansellhealthcare.com

www.ansell.com

LEADERS IN HEALTHCARE SAFETY THROUGH COLLABORATION AND INNOVATION

SANDEL healthcare safety devices provide an array of safety solutions to protect patients and healthcare workers. These innovative safety products were developed through collaboration with nurses and surgical teams by listening to their needs and suggestions.

To submit your ideas, please visit: www.sandelmedical.com/ideas

References:

1. Centers for Disease Control and Prevention. National Institute for Occupational Safety and Health (NIOSH). NIOSH Alert: Preventing Needlestick Injuries in Health Care Settings. Published November 1999. NIOSH publication 2000-108. <http://www.cdc.gov/niosh/docs/2000-108/pdfs/2000-108.pdf>. Accessed February 22, 2012. **2.** Centers for Disease Control and Prevention. National Institute for Occupational Safety and Health (NIOSH). Slip, Trip, and Fall Prevention for Healthcare Workers. Published December 2010. NIOSH publication 2011-123. www.cdc.gov/niosh/docs/2011-123/pdfs/2011-123.pdf. Accessed February 23, 2012. **3.** US Department of Labor. Occupational Safety & Health Administration. Foundry ergo partnership (FEP). https://www.osha.gov/dcsp/partnerships/regional/region5/722_20030314_fep.html. Accessed March 8, 2012. **4.** US Food & Drug Administration. Med error reports to FDA show a mixed bag. <http://www.fda.gov/downloads/Drugs/DrugSafety/MedicationErrors/ucm115775.pdf>. Published October 1, 2001. Accessed March 8, 2012. **5.** US Food & Drug Administration. Medication error reports. <http://www.fda.gov/Drugs/DrugSafety/MedicationErrors/ucm080629.htm>. Accessed March 8, 2012. **6.** Wagar EA, Stankovic AK, Raab S, Nakhleh RE, Walsh MK. Specimen labeling errors. A Q-probes analysis of 147 clinical laboratories. *Arch Pathol Lab Med.* 2008; 132:1617-1622. **7.** World Health Organization. *WHO Guidelines for Safe Surgery 2009: Safe Surgery Saves Lives.* Geneva, Switzerland: World Health Organization; 2009.

* and ™ are trademarks owned by Ansell Limited or one of its affiliates.
© 2012 Ansell Limited. All rights reserved. 130326 Rev. A

